

The Asterisk Cell Knee Protection System is without question, the most advanced-designed knee system available. Patent pending revolutionary design concepts in it's key features are apparent at first site and confirmed by it's first use. Born from riders with orthopedic bracing backgrounds, the off-the-shelf Cell becomes a custom-fit device though the riders own simple adjustments. Simply, the Cell is unchallenged as the best off-the-shelf knee protection available, answering the need for superior protection at a reasonable price.

Size	(in)	(cm)
Small	13" - 13 7/8"	33.02 - 35.24
Medium	14" - 14 7/8"	35.56 - 37.78
Large	15" - 15 7/8"	38.1 - 40.32
X-Large	16" - 16 7/8"	40.64 - 42.86
XX-Large	17" - 18 1/4"	43.18 - 46.36

How to size the Asterisk Cell Knee Brace:

The circumference chart is a guideline only. This chart is used to help indicate the appropriate size. The most important indicator is the Mid Kneecap measurement. This measurement should be the basis for the size. To ensure a proper fit, circumference measurements should be taken directly on the skin with the leg straightened, using a soft tape measure. (If a soft tape measure is not available, a piece of string or yarn may be used to calibrate the leg's diameter, then measured against a hard tape measure).



Locate the middle of the kneecap and measure around the leg. Measurements should be taken directly on the skin with a straight leg.